

Tip of the Month: Troubleshoot MaxiMove Dependent Lift




Problem: MaxiMove does not turn on when GREEN power button is pressed. **Check the following:**

1. Battery life indicator on the hand control; replace battery if low or not charged.



Full = charged
Low = needs charging



2. Spare battery is charging properly; if not, enter Clinical Engineering work order for charger.

Amber light on charger =
Battery charging



Green light on charger =
Battery fully charged



3. If battery is replaced and MaxiMove still does not work:

- Enter Clinical Engineering work order for MaxiMove; print work order and tape to lift
- Label malfunctioning lift, battery and/or battery charger
- Communicate to manager and colleagues