

Tip of the Month: Troubleshoot SaraPlus Sit-Stand Lift



Problem: SaraPlus does not turn on when the GREEN power button is pressed. **Check:**

1. Emergency Stop on the side of SaraPlus; if pressed in, twist and pull out to deactivate



Emergency Stop Deactivated:
Lift ready to use



Emergency Stop Activated:
Red top pushed in

2. Battery charge on Battery Indicator on the right side of lift; replace battery if low or not charged



Battery Indicator
Green = charged
Yellow = low charge/need to charge

Problem: Battery is replaced, but SaraPlus still does not turn on. **Check:**

The battery charger; if charging properly, you will see:



Battery Charger:
Amber light =
charging



Battery Charger:
Green light = battery
charged

Problem: All previous steps have been taken and SaraPlus still does not turn on. **Action:**

- Enter Clinical Engineering work order for SaraPlus lift
- Print work order and tape to device
- Label malfunctioning lift, battery and/or battery charger
- Communicate to manager and colleagues