

## CHECKLIST FOR INSPECTION OF EYEWASH OFF SITE CLINICS

1. Ensure that the path to the eyewash/drench hose is not obstructed.
2. Verify that nozzle caps are in place to prevent contamination and that the nozzles, nozzle caps, and bowl/sink are clean and sanitary.
3. Actuate valve to full open position. Water must flow within 1 second.
4. Verify that nozzle caps come off when the eyewash or drench hose is activated.
5. Verify that water continues to flow until manually turned off and can be used without requiring the use of the operator's hands.
6. Look at the flow pattern. It should provide a gentle non-injurious flow. If a dual-stream eyewash, both streams should rise to equal height in a pattern that will flush both eyes simultaneously.
7. Continue to flush until water is clear.
8. Report problems to your building maintenance provider.
9. **Annual Check:** Temperature \_\_\_\_\_F Flow rate checked Met Not Met (repair date \_\_\_\_\_) Initials \_\_\_\_\_

Month	Date and initial to show above items have been checked					Note Problems (reported/resolved)
January	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
February	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
March	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
April	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
May	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
June	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
July	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
August	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
September	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
October	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
November	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	
December	Week 1	Week 2	Week 3	Week 4	Week 5 or NA	