



From stretcher:

1. Raise the head of stretcher up between 45 to 90 degrees.
2. Place a blue liner behind the seated patient, then slide the sling down between the blue liner and patient's back.
3. Position Maximove from the side of the stretcher, open the base legs, if needed.

From wheelchair:

1. Remove arm rests and foot rest if able.
2. Place a blue liner behind the seated patient, then slide the sling down between the blue liner and patient's back.
3. Open Maximove's base legs and straddle the wheelchair at an angle,
4. Turn the hanger bar and align it to the patient.