

# WORKING SAFELY WITH ANESTHETIC GASES IN THE LABORATORY

## TIPS FOR SAFE ANESTHETIC GAS USE

### EXAMPLES OF VETERINARY ANESTHESIA EQUIPMENT

- Mouse breathing circuit from Harvard Apparatus



- VetEquip VaporGuard Charcoal canister



### OCCUPATIONAL HYGIENE AND SAFETY DIVISION

Box 3387 DUMC  
Durham, NC 27710  
Phone: 919-684-5996  
Fax: 919-681-5916

### EQUIPMENT AND WORK PRACTICES

- Use a vaporizer (or anesthesia machine) and breathing circuit/nose cone designed for anesthesia for the animal species you are working with. (Do NOT attempt to “improvise” anesthesia delivery with equipment not created for that purpose.)
- If using an induction chamber, flush the chamber before opening (with air or O<sub>2</sub>) and/or minimize the time the chamber is open.

### SCAVENGING WASTE GASES

- Work in a fume hood if possible or near a dedicated local exhaust “snorkel” positioned 4–6” from the induction chamber or nose cone.
- Medical vacuum systems may be used to exhaust waste anesthetic gases from an anesthesia vaporizer circuit. If your building has a lab vacuum system, it is designed for air only—check with OESO for approval to use a standard lab vacuum system for scavenging.
- If using an **isoflurane vaporizer**, you can work on the benchtop if you have a charcoal canister for scavenging. Follow the manufacturer’s instructions for setting up the canister and tracking canister life.

### FOR MORE INFORMATION

- See the Duke Animal Care and Use Program’s policy on [Safe Use of Anesthetic Gases](#).
- If you are concerned about employee exposure to anesthetic gases, contact OESO at 919-684-5996.
- For guidelines on animal safety during anesthesia, see the website of the [Duke Animal Care and Use Program](#) or contact DLAR or OAWA.