

## Guide to Estimating Amputee Sling Sizes for Patient Height and Weight

**Note:** When choosing sling size, consider a patient's height, weight, width, and girth.  
This is **ONLY A GUIDE** as body types vary.

| Height / Weight | < 90 lbs/<br><41 kg | 90-139 lbs/<br>41 – 63 kg | 140-219 lbs/<br>64 – 99 kg | 220-500 lbs/<br>100 – 226 kg | 501 – 600 lbs/<br>227 – 272 kg                           | 601 - 800 lbs/<br>273 – 363 kg  |
|-----------------|---------------------|---------------------------|----------------------------|------------------------------|--|---|
| 3'6" - 4'5"     | M/L                 | M/L                       | M/L                        | XL/XXL                       | XL/XXL<br>But <u>ONLY</u> for<br>use with<br>MaxiSky 600 | XL/XXL<br>But <u>ONLY</u> for<br>use with<br>MaxiSky 1000 &<br>SizeWise Lifts |
| 4'6" - 5'4"     | M/L                 | M/L                       | M/L                        | XL/XXL                       |  |   |
| 5'5" - 5'11"    | M/L                 | M/L                       | M/L                        | XL/XXL                       |  |   |
| 6'0" - 6'4"     | M/L                 | XL/XXL                    | XL/XXL                     | XL/XXL                       |  |   |
| > 6'5"          | XL/XXL              | XL/XXL                    | XL/XXL                     | XL/XXL                       |  |   |


**MaxiMove/Opera/Tempo**  
Sling with CLIP attachment



**MaxiSky**  
Sling with LOOP attachment



*Amputee  
Option*



Leg wing/straps slide  
into the amputee strap  
on the opposite leg



Updated 01-20