Duke Health- Bariatric Equipment Rentals

Equipment order and Discharge process:

- 1. Follow Bariatric Plan of Care
- 2. Select a bed and mattress type based on the patient's weight, height, and girth.
- 3. Select other ancillary items as needed per BMAT level.
- 4. To order and return equipment at discharge, call Agiliti (Sizewise) @ 1-800-814-9389
 - DUH: 1) Online "Request of Beds and Equipment Distribution" or 2) call Agiliti Health (Sizewise)
 - DRH: call Agiliti Health (Sizewise)
 - DRAH: call Patient Care Equipment , Cheronte Ross @ 919-862-5850
 - Ambulatory Clinics: call Agiliti Health (Sizewise)

NOTE: Return any Agiliti (Sizewise) equipment immediate after patient is discharged to avoid unnecessary charges to the unit or department.

Available Equipment from Agiliti (Sizewise) and Descriptions:



Bed Platform: BariRehab 3

- 1000 lb. capacity
- Width: 39" (for < 500 lbs/227 kg) or 48" (for > 500 lbs/227 kg)
- Length: 80" (for under 6'2") or 86" (for over 6' 2")
- Lowest bed height: 15"
- Power-drive
- Scale
- Bed exit alarm
- Trapeze option

Mattress Option #1: Big Turn 2 (Typically used in ICU)

Low air loss

Optional features: pulsation therapy, or full body lateral rotational therapy with percussion and vibration

Mattress Option # 2: Mighty Air (Typically used in step-down)

• Low air loss & alternation therapy



Bed Platform: Expandable Evolution

- 850 lbs/364 kg capacity
- Width: 39" (for < 500 lbs/227 kg) or 48" (for > 500 lbs/227 kg)
- Length: 82"
- Lowest bed height: 9"
- Scale
- Bed exit alarm
- Trapeze option

Mattress Options:

Low air loss

Optional features: pulsation therapy, or full body lateral rotational therapy with percussion and vibration

• Foam (ONLY at DRH & DRAH for bariatric by-pass surgical patients)



"Shuttle" Chair/Stretcher

Option 1: Shuttle A Advanced (650 lb capacity) — Seat Width: 25.5"
Option 2: Shuttle B Advanced (1000 lb capacity) - Seat Width: 31"

Features:

- Height in stretcher mode for both models: 27.5" 35.5"
- Power-operated
- Chair and stretcher positions
- Seat assist to lift patient to stand
- Independently elevating leg rests

Updated: Aug 2023 Page 1 of 2

	 Trapeze 850 lbs/386 kg capacity Available for Bari Rehab3 and Expandable Evolution
	 Bari Lift and Transfer, Bedside 750 lbs/341 kg or 1000 lbs. capacity Power-operated Scale option
	Sit-to-stand lift • 800 lbs • Power-operated • Lifts to stand and hold in standing • Transfers between bed, bedside commode, wheelchair, and recliner
	Bedside Recliner Too lbs. capacity Seat Width: 29" Height: 21"
	 Bari Drop Arm Commode 600 lbs/273 kg, Seat width 24" 850 lbs/386 kg, Seat width 28" *Height: 18" – 22" (adjustable legs for both sizes above) 1000 lbs/455 kg, Seat width 36", non-height adjustable
ARA	 Bari-Shower Commode 750 lb. capacity Seat Width: 26" or 30"; Height: 22" Removable back, arms, and swing-away footrests
	 Bari Walker 750 lb. capacity Width: 27" Handle Height: 31" - 39 "
	 Bari Wheelchair 850 lb. capacity Seat Width: 24", 26", 28" (commonly rented) Other Width available: 30", 32", 34', 36" Removable Armrests Elevating Leg Rests (Swing Away or Removable)

Updated: Aug 2023 Page 2 of 2